CONTENT

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Tuloc®

Tumor Localizer

The adjustable TULOC® Tumor Localizer kit from Somatex® is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate intraoperative localization of the findings by the surgeon. The development of this product complies with the latest engineering knowledge featuring the highest level of comfort, quality and reliability for both the physicians and the patients.

The TULOC® Tumor Localizer consists of a hypodermic needle and a preloaded marking wire with distal double arches (Fig. 1).

The following steps must be adhered to prior to the application:

1. Check that the package has not been opened or damaged before use; check the sterility date.
2. Open the package.

Note:

The shaft and the double arches of the marking wire are situated within the hypodermic needle in the delivered condition. Before using the marking wire, ensure that both of the first dark annular markings behind the needle handle are visible (Fig. 3, (3)). In this position, the arches of the wire are positioned completely within the needle.

Description of application:

3. Local disinfection and administration of local anesthetic in the insertion area.
4. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
5. Needle insertion insert the needle under sonographic or mammographic guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
6. After achieving the target point, you can place a rolled wire into the tumor area (Fig. 2, (2)). To achieve this, you should hold the needle position and introduce the wire into the cannula in the direction from the proximal to the distal point. Then the rolled wire should be completely deployed, as soon as the marking reaches the needle insertion (Fig. 2, (1)).
Note: The wire may be pushed into the cannula only until the last visible marking on the wire (1) has been inserted into the cannula attachment!
7. Prior to the removal of the needle, a proper positioning of the rolled wire can be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the needle position, the wire can be deployed again according to the description in point 5. and 6. for preoperative marking.
8. Remove the needle from the mamma carefully.
9. Relax the compression of the breast, if necessary. Push the opened fixation part in the compressed state of the breast above the wire up to the surface of the skin (**Caution:** some space must be left between the skin and the fixation part). The wire must lie between the brackets (Fig.4, (4)) of the fixation part. The wire is fixed through the closure being locked in place (Fig. 4, (5)). This measure prevents the risk of change in position of the wire along the distal end.
10. Treat the wound.
11. Dispose of the puncture needle used.

Contra-indications:

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply for the use of the retractable Somatex®TULOC® Tumor Localizer.

Particularly for Applications in Magnetic Resonance Imaging

In the event that electroconductive materials are inserted into the body at 1.5 T* >10 cm at 3 T >5 cm, there is a possibility that the tip of the conductor or the cannula, respectively, will heat to a temperature in excess of 70 °C. The occurrence of this effect can be minimized by using low-energy sequences.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilisation or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilisation processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

Storage instructions:

Store in a dark, dry place (humidity 30 – 65%) at a temperature between 41 and 86°F (5 – 30°C).

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Somatex Medical Technologies GmbH