
TUMARK[®] Professional

zur perkutanen Tumormarkierung
for percutaneous tumor marking

SOMATEX[®]

MEDICAL TECHNOLOGIES GMBH

MADE IN GERMANY

STERILE EO

CE 0123

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INSTRUCTION OF USE

Area of Application:

The Tumark[®] serves for percutaneous marking in soft tissue, in particular in the pectoral tissue. Among other things, the markers are implanted in lesions prior to or during chemotherapy, the location of a biopsy sampling point or a removed tumor, as well as for enhanced orientation for planning radiation treatment.

Note

This product has been explicitly designed and manufactured for single use only. The product may not be re-used after a single application. The quality of materials, coatings and adhesive joints may deteriorate. A secure application is no longer guaranteed. The product is not designed for cleaning and sterilization processes after a single use. The sterility of reprocessed single-use is therefore not guaranteed. The risk of unintentional injuries and infections, especially cross-infection with a patient and medical personnel increases unmeasurably. The company Somatex will not be held responsible in any way for the consequences in case of resterilisation. The nickel-titanium (Nitinol) marker (3) preloaded into the distal end of the cannula (1) is not absorbable and only intended for one-time use. In the event of the reuse of even one of the components of the product, Somatex will deny any responsibility. The marker itself is always MR-compatible, the puncturing tool according to the order number.

Application Description:

1. Prior to use, make sure that the package has not been opened and/or damaged; check the sterilization date,

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* Tesla

2. Disinfection of the puncturing point,
3. In an imaging procedure such as ultrasound, Mammogram, or when using a MR-compatible set, in magnetic resonance imaging as well, puncture the target area with the cannula (1).
4. Place the marker (3) by pushing the slide (2) as far as it will go,
5. Document the location of the marker (3),
6. Remove the cannula (1),
7. Treat the wound,
8. After the treatment, please ensure the appropriate disposal of the cannula (1) in the proper cannula container.

Contraindications:

For the application of the cannulae, all the contraindications for the respective area of application apply as they are understood to be known for cannulae of this type in accordance with the rules of medicine.

Particularly for Applications in Magnetic Resonance Imaging

In the event that electroconductive materials are inserted into the body at 1.5 T* >10 cm at 3 T >5 cm, there is a possibility that the tip of the conductor or the cannula, respectively, will heat to a temperature in excess of 70 °C. The occurrence of this effect can be minimized by using low-energy sequences.

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